



for better life

A blurred photograph of a young child lying down, possibly in a hospital bed, wearing a blue medical device around their neck. The child's face is in focus but the background is out of focus.

Respiratory & Anesthesia



Respiratoire & Anesthésie



OUR MISSION:

Since the moment of establish, Henso are committed to offer one-stop solution for medical supplies, to reduce medical expense by supply reliable products with affordable cost. We are efforts for human's better life.

GOOD MANUFACTURING PRACTICES:

Our Manufacturing plant conforms to GMP, a WHO proprietary benchmark. This standard is evident in the physical ambience of our plant, the air and water quality, manufacturing technology, packaging and the sterilization process. Ambience: Special micro and HEPA filter at all air handling units ensure minimal contamination in all manufacturing areas. Micro organism and particle concentration levels are carefully monitored by the micro-biological lab to keep these levels within safe limits. Employees undergo periodic medical checkup to maintain high standards of health and hygiene. All clean rooms are provided with a three step change room and air showers at the entrance to ensure a dust-free environment.

PURSUIT ZERO DEFECT QUALITY:

Manufacture zero-defect products is our goal of quality. We have built a series online and manual quality control measures into each stage of raw material procurement, manufacturing and packing cycle. Our Quality Control standards conform to ISO 13485:2012 and European Medical Device Directive 93/42/EEC Annex V. We also registered with US FDA. As a result, our products have been authorized to lable prestigious CE mark.

PACKING AND STERILIZATION:

Our products are packed by most modern packaging machines, to ensure minimal risk of contamination or damage during transportation. We use Ethylene Oxide to sterilize the product, as per ISO 11135/EN 550 standards, each batch were release after strict inspection.



Venturi Mask with 2 diluters
Masque Venturi avec 2 diluteurs

HSA0101

- Made from medical grade PVC
- Adjustable nose clip
- Anti-crush tube with 2m length
- Green or transparent mask & tube
- Oxygen Concentration:
Green diluter (24%, 26%, 28% and 30%)
White diluter (35%, 40% and 50%)
- Size: S, M, L, XL

Venturi Mask with 6 diluters
Masque Venturi avec 6 diluteurs

HSA0102

- Made of medical grade PVC
- Adjustable nose clip
- Anti-crush tube with 2m length
- Green or transparent mask & tube
- Oxygen Concentration: blue 24%, yellow 28%, white 31%, green 35%, pink 40%, orange 50%
- Size: S, M, L, XL



Aerosol Mask
Masque Aérosol

HSA0103

- Made of medical grade PVC
- Adjustable nose clip
- Available in Green or transparent color
- DEHP free and 100% latex free available
- Size: S (pediatric standard), M (pediatric elongated), L (adult standard), XL (adult elongated)



Nebulizer Mask

Masque Nébuliseur

HSA0104

- Made of medical grade PVC
- Adjustable nose clip
- Anti-crush tube with 2m length
- Green or transparent mask & tube
- Optional nebulizer with three types
- Non-sterile or Sterile EO, individual pack
- Size: S, M, L, XL

Nebulizer Kit

Kit nébuliseur

HSA0105

- Made of medical grade PVC
- Anti-crush tube with 2m length
- Green or transparent color
- Optional nebulizer with three types A/B/C
- Nebulizer capacity: 6ml, 8ml, 10ml
- Non-sterile or Sterile EO, individual pack




Oxygen Mask

Masques à Oxygène

HSA0106

- Designed for convenient oxygen therapy
- Made of medical grade PVC, non-toxic, latex free
- Medium concentration, adjustable nose clip
- Available in green or transparent color
- Non-sterile or Sterile EO, individual pack
- Size: S (pediatric standard), M (pediatric elongated), L (adult standard), XL (adult elongated)





High Concentration Oxygen Mask | HSA0107

Masque à oxygène haute concentration

- Designed for high concentration oxygen therapy
- Made of medical grade PVC, non-toxic, latex free
- High concentration, adjustable nose clip
- Anti-crush tube with 2m length
- Green or transparent mask and tube
- Non-sterile or Sterile EO, individual pack
- Size: S, M, L, XL

Anesthesia Mask | HSA0108

Masques d'anesthésie

- Made of medical grade PVC
- Air cushion guarantee comfortable and safety seal
- Color coded valve for easy size identification
- Clear and transparent mask
- Non-sterile, individual packed in a poly bag

Size	0#	1#	2#	3#	4#	5#
Patient	Neonate	Infant	Pediatric	Adult Small	Adult Medium	Adult Large



CPR Mask | HSA0109

Masque CPR

- Made of medical grade PVC, non-toxic, latex free
- Prevents direct contact during artificial respiration
- Transparent dome allows rescuer to visually check for patient's lip colour and vomitus
- Plastic hard box in red color
- Non-sterile, individual packed in a poly bag
- Size: Pediatric, Adult



CPR Life Key Life Key Masque

HSA0110

- Designed for emergency first aid
- Prevents direct contact during rescuer and victim
- Nylon bag + PVC mask
- Surrounded by a transparent face sheet that covers the whole face as an additional protection
- PVC Mask size: 27.5 x 20cm

Disposable PVC Laryngeal Mask Airway Masque laryngé à usage unique PVC

HSA0111

- Made of medical grade PVC
- Color coded for easy size identification
- Sterile EO, for single use only

Size (#)	1.0	1.5	2.0	2.5	3.0	4.0	5.0
Inflation Volume (ml)	4	7	10	14	20	30	40
Patient Weight (KG)	0-5	5-10	10-20	20-30	30-50	50-70	70-100



Disposable Silicone Laryngeal Mask Airway Masque laryngé à usage unique Silicone

HSA0112

- Made of medical grade Silicone
- Silicone soft cuff provides good sealing
- Cuff color optional in Flesh, Blue or Clear

Size (#)	1.0	1.5	2.0	2.5	3.0	4.0	5.0
Inflation Volume (ml)	4	7	10	14	20	30	40
Patient Weight (KG)	0-5	5-10	10-20	20-30	30-50	50-70	70-100





Reinforced Silicone Laryngeal Mask Airway

Masque laryngé renforcé à usage unique Silicone

HSA0113

- Made of medical grade PVC
- Reinforced by stainless steel spring
- Sterile EO, for single use only

Size (#)	1.0	1.5	2.0	2.5	3.0	4.0	5.0
Inflation Volume (ml)	4	7	10	14	20	30	40
Patient Weight (KG)	0-5	5-10	10-20	20-30	30-50	50-70	70-100

Reusable Silicone Laryngeal Mask Airway

Masque laryngé réutilisable Silicone

HSA0114

- Made of medical grade Silicone
- Autoclave sterilization only, resistance at 134°C
- Reusable up to 40 times

Size (#)	1.0	1.5	2.0	2.5	3.0	4.0	5.0
Inflation Volume (ml)	4	7	10	14	20	30	40
Patient Weight (KG)	0-5	5-10	10-20	20-30	30-50	50-70	70-100



Reusable Reinforced Silicone Laryngeal Mask Airway

Masque laryngé renforcé réutilisable Silicone

HSA0115

- Made of medical grade Silicone
- Reinforced by stainless steel spring
- Reusable up to 40 times with record card

Size (#)	1.0	1.5	2.0	2.5	3.0	4.0	5.0
Inflation Volume (ml)	4	7	10	14	20	30	40
Patient Weight (KG)	0-5	5-10	10-20	20-30	30-50	50-70	70-100



Berman Oropharyngeal Airway | HSA0116

Canules de berman

- Made of non-toxic polyethylene, latex free
- Berman type, smooth surface
- Color coded for easy size identification
- Sterile EO, individual pack in poly bag

Size	000	00	0	1	2	3	4	5	6
Length (mm)	40	50	60	70	80	90	100	110	120
Color	Pink	Blue	Black	White	Green	Yellow	Red	Light Blue	Orange

Guedel Oropharyngeal Airway | HSA0117

Canules de guedel

- Made of non-toxic polyethylene, latex free
- Guedel type, smooth surface
- Color coded for easy size identification
- Sterile EO, individual pack in poly bag

Size	000	00	0	1	2	3	4	5	6
Length (mm)	40	50	60	70	80	90	100	110	120
Color	Pink	Blue	Black	White	Green	Yellow	Red	Light Blue	Orange



Nasopharyngeal airway | HSA0118

Canules nasopharyngé

- Special designed for nasal intubation
- Made of medical grade PVC, non-toxic, latex free
- Atraumatic soft rounded beveled tip
- Blue, Green or white semi-transparent color
- Sterile EO for single use, individual pack
- Size (mm): 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0



Disposable PVC Endotracheal Tube

Sonde endotrachéale à usage unique PVC

HSA0119


- Made of soft and clear PVC, non-toxic
- Size and mark clearly printed on tube
- X-ray opaque line and Murphy eye design
- Available in with cuff or without cuff
- Sterile EO, for single use only
- Size (mm): 2.0, 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10.0

PVC Reinforced Endotracheal Tube

Sonde endotrachéale armée à usage unique PVC

HSA0120

- Made of medical grade PVC
- Reinforced by stainless steel spring
- High volume and low pressure cuff
- X-ray opaque line and Murphy eye design
- Size and mark clearly printed on tube
- Sterile EO, for single use only
- Size (mm): 2.0, 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10.0




Silicone Reinforced Endotracheal Tube

Sonde endotrachéale armée à usage unique Silicone

HSA0121

- Made of silicone, soft and clear
- Reinforced by stainless steel spring
- High volume and low pressure silicone cuff
- X-ray opaque line and Murphy eye design
- Size and mark clearly printed on tube
- Sterile EO, for single use only
- Size (mm): 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0



Reusable Silicone Reinforced Endotracheal Tube

Sonde endotrachéale armée réutilisable Silicone

HSA0122

- Made of silicone, soft and clear
- Reinforced tube to prevent crush and kinking
- High volume and low pressure silicone cuff
- X-ray opaque line and Murphy eye design
- Size and mark clearly printed on tube
- Reusable up to 10 times
- Size (mm): 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0

Oral Preformed Endotracheal Tube

Sonde endotrachéale préformes orale avec ballonnet

HSA0123

- Special designed for oral intubation
- Made of non-toxic PVC, latex free
- Oral preformed tube with X-ray opaque line
- With Murphy eye design
- Sterile EO, individual pack for single use only
- Size (mm): 2.0, 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10.0





Nasal Preformed Endotracheal Tube

Sonde endotrachéale préformes nasale avec ballonnet

HSA0124

- Special designed for nasal intubation
- Made of non-toxic PVC, latex free
- Nasal preformed tube with X-ray opaque line
- With Murphy eye design
- Sterile EO, individual pack for single use only
- Size (mm): 2.0, 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10.0






Endotracheal Tube (Suction Type) | **HSA0125**
 Sonde endotrachéale (Type d'aspiration)

- With extra suction tube or suction catheter
- Made of soft and clear PVC, non-toxic
- Size and mark clearly printed on tube
- X-ray opaque line and Murphy eye design
- Available in with cuff or without cuff
- Sterile EO, individual pack for single use only
- Size (mm): 2.0, 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10.0

Endotracheal Tube with Mucus Extractor
 Sonde endotrachéale (Type Extracteurs de mucus)

HSA0126

- With extra Mucus Extractor
- Size and mark clearly printed on tube
- X-ray opaque line and Murphy eye design
- Available in with cuff or without cuff
- Sterile EO, individual pack for single use only
- Size (mm): 2.0, 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10.0



Tracheostomy Tube | **HSA0127**
 Canule de trachéotomie

- With soft-double wings and neck tie
- Less invasive atraumatic soft rounded tip
- X-ray line help precisely location
- Pre-curved shape easy to insert, reduce trauma
- Available in With cuff or without cuff
- Sterile EO, individual pack for single use
- Size(mm): 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10.0



Tracheostomy Tube (Suction Type)

Canule de trachéotomie (Type d'aspiration)

HSA0128


- With extra suction tube or catheter
- Less invasive atraumatic soft rounded tip
- X-ray line help precisely location
- Pre-curved shape easy to insert, reduce trauma
- With or without cuff
- Sterile EO, individual pack for single use
- Size(mm): 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10.0

Tracheostomy Tube with Mucus Extractor

Canule de trachéotomie (Type Extracteurs de mucus)

HSA0129

- With Mucus Extractor & Suction Catheter
- Less invasive atraumatic soft rounded tip
- X-ray line help precisely location
- Pre-curved shape easy to insert, reduce trauma
- Sterile EO, individual pack for single use
- Size(mm): 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10.0



Suction Connecting Tube

Tube d'aspiration de connexion

HSA0130

- Made from clear medical grade PVC
- Anti-kinking tube to avoid blocking
- Connector tip available in four types
- Tube Diameter: 1/4" or 3/16"
- Tube Length: 1.8m, 2.5m, 3.6m or customized
- Connector types: A, B, C, D
- Sterile EO, individual pack for single use





Suction Connecting Tube with Yankauer Handle
Tube d'aspiration de connexion Yankauer Poignée

HSA0131

- Connect with Yankauer Handle
- Made of clear medical grade PVC
- Anti-kinking tube to avoid blocking
- Connector tip available in four types
- Tube Diameter: 1/4" or 3/16"
- Tube Length: 1.8m, 2.5m, 3.6m or customized
- Sterile EO, individual pack for single use

Yankauer Handle

Canule d'aspiration Yankauer

HSA0132

- Made of medical grade k-resin, non-toxic
- Used together with Suction Connecting tube
- Crown Tip or Standard (flat) Tip
- With control vent or without control vent
- Flash surface, 100% transparent
- Edge polished and trimmed, trauma free
- Sterile EO, individual pack for single use



Anesthesia Breathing Circuit

Circuit d'anesthésie de respiration

HSA0133

- Made of transparent medical grade EVA
- Standard connector: Pediatric 15mm, Adult 22mm
- Kinking resistant tube to avoid blocking
- Match all kinds of anesthesia machine
- Length: 1.0m, 1.2m, 1.5m, 1.6m, 1.8m, 2.0m or customized



Anesthesia Breathing Circuit (Expandable)
Circuit d'anesthésie de respiration extensible

HSA0134

- Expandable hoses are easy for transportation
- Made of transparent medical grade EVA
- Standard connector: Pediatric 15mm, Adult 22mm
- Kinking resistant tube to avoid blocking
- Match all kinds of anesthesia machine
- Length: 1.0m, 1.2m, 1.5m, 1.6m, 1.8m, 2.0m

Anesthesia Breathing Circuit Set (A)
Circuit d'anesthésie de respiration ensemble

HSA0135

- Made of medical grade plastic material
- Corrugated hoses + Breathing bag + Anesthesia Mask + HME Filter
- Standard connector: Pediatric 15mm, Adult 22mm
- Kinking resistant tube to avoid blocking
- Length: 1.0m, 1.2m, 1.5m, 1.6m, 1.8m, 2.0m



Anesthesia Breathing Circuit Set (B)
Circuit d'anesthésie de respiration ensemble

HSA0136

- Made of medical grade plastic material
- Corrugated hoses + Breathing bag + Anesthesia Mask + Bacterial Filter
- Standard connector Pediatric 15mm, Adult 22mm
- Kinking resistant tube to avoid blocking
- Length: 1.0m, 1.2m, 1.5m, 1.6m, 1.8m, 2.0m

Quality Assurance

EC Certificate

Production Quality Assurance System
Directive 93/42/EEC on Medical Devices (MDD), Annex V
(Devices in Class A, B, or C)
No. Q2 14 09 89326 002

Manufacturer: Hensol Medical (Hangzhou) Co., Ltd.
Room 1202, 4A Cambridge Court
No. 100 Shuang Road
31000 Hangzhou
PEOPLE'S REPUBLIC OF CHINA

EC-Representative: Weikang Ltd
Suite 8, 28 Halsey Street
LONDON
W1D 2JH
UNITED KINGDOM

Product Category(ies): For detail information see attachment

The Certification Body of TÜV SÜD Product Service GmbH declares that the aforementioned manufacturer has implemented a quality assurance system for production and final inspection of the respective devices in accordance with MDD Annex V. This quality assurance system conforms to the requirements of this Directive and is subject to periodic surveillance. For expiring of class B and C devices an additional Annex IX certificate is mandatory. See also notes overleaf.

Report No.: 041108001

Valid from: 2015-01-13
Valid until: 2020-01-12

Date: 2015-01-14

TÜV SÜD Product Service GmbH is notified body with identification no. 0123
Page 1 of 3

TÜV SÜD Product Service GmbH - Zertifizierungsstelle - Referenz-Nr.: 8028 München - Germany



EC Certificate

Production Quality Assurance System
Directive 93/42/EEC on Medical Devices (MDD), Annex V
(Devices in Class A, B, or C)
No. Q25 14 09 89326 003

Manufacturer: Hensol Medical (Hangzhou) Co., Ltd.
Room 1202, 4A Cambridge Court
No. 100 Shuang Road
31000 Hangzhou
PEOPLE'S REPUBLIC OF CHINA

EC-Representative: Weikang Ltd
Suite 8, 28 Halsey Street
LONDON
W1D 2JH
UNITED KINGDOM

Product Category(ies): Sterile Urine Drainage Bag, Sterile Vaginal Speculum, Disposable Umbilical Cord Clamps, Disposable I.D. Brackets, Disposable Dressing Kits, Disposable Spermocidal Suits, Gauze Swab (Sponges), Gauze Bandages, Absorbent Cotton Balls, Absorbent Cotton Balls, First Aid Kits, Bandages, Sterile Wound Dressing, Sterile Wound Dressing, Surgical Gowns, Face Masks, Shoe Covers, Cotton Swab (Cotton Tipped Applicators), Tongue Depressor, Sterile Feeding Tubes, Disposable Stomach Tubes, Disposable Stomach Tubes, Sterile Latex Examination Gloves

The Certification Body of TÜV SÜD Product Service GmbH declares that the aforementioned manufacturer has implemented a quality assurance system for production in accordance with MDD Annex V. This quality assurance system covers those aspects of manufacture concerned with ensuring and maintaining the conformity of the respective devices / device categories and conforms to the requirements of this Directive. It is subject to periodic surveillance. See also notes overleaf.

Report No.: 041108001

Valid from: 2015-01-13
Valid until: 2020-01-12

Date: 2015-01-14

TÜV SÜD Product Service GmbH is notified body with identification no. 0123
Page 1 of 3

TÜV SÜD Product Service GmbH - Zertifizierungsstelle - Referenz-Nr.: 8028 München - Germany



EC Certificate

Production Quality Assurance System
Directive 93/42/EEC on Medical Devices (MDD), Annex V
(Devices in Class A, B, or C)
No. Q2M 14 09 89326 004

Manufacturer: Hensol Medical (Hangzhou) Co., Ltd.
Room 1202, 4A Cambridge Court
No. 100 Shuang Road
31000 Hangzhou
PEOPLE'S REPUBLIC OF CHINA

EC-Representative: Weikang Ltd
Suite 8, 28 Halsey Street
LONDON
W1D 2JH
UNITED KINGDOM

Product Category(ies): Aneroid Sphygmomanometers

The Certification Body of TÜV SÜD Product Service GmbH declares that the aforementioned manufacturer has implemented a quality assurance system for the manufacture in accordance with MDD Annex V. This quality assurance system covers those aspects of manufacture concerned with the manufacturing requirements of the respective devices / device categories and conforms to the requirements of this Directive. It is subject to periodic surveillance. See also notes overleaf.

Report No.: 041108001

Valid from: 2015-01-13
Valid until: 2020-01-12

Date: 2015-01-14

TÜV SÜD Product Service GmbH is notified body with identification no. 0123
Page 1 of 3

TÜV SÜD Product Service GmbH - Zertifizierungsstelle - Referenz-Nr.: 8028 München - Germany



EC Certificate

Production Quality Assurance System
Directive 93/42/EEC on Medical Devices (MDD), Annex V
(Devices in Class A, B, or C)
No. Q2 14 09 89326 002

Manufacturer: Hensol Medical (Hangzhou) Co., Ltd.
Room 1202, 4A Cambridge Court
No. 100 Shuang Road
31000 Hangzhou
PEOPLE'S REPUBLIC OF CHINA

EC-Representative: Weikang Ltd
Suite 8, 28 Halsey Street
LONDON
W1D 2JH
UNITED KINGDOM

Product Category(ies): For detail information see attachment

The Certification Body of TÜV SÜD Product Service GmbH declares that the aforementioned manufacturer has implemented a quality assurance system for production and final inspection of the respective devices in accordance with MDD Annex V. This quality assurance system conforms to the requirements of this Directive and is subject to periodic surveillance. For expiring of class B and C devices an additional Annex IX certificate is mandatory. See also notes overleaf.

Report No.: 041108001

Valid from: 2015-01-13
Valid until: 2020-01-12

Date: 2015-01-14

TÜV SÜD Product Service GmbH is notified body with identification no. 0123
Page 1 of 3

TÜV SÜD Product Service GmbH - Zertifizierungsstelle - Referenz-Nr.: 8028 München - Germany



Attachment for Certificate No. Q25 14 09 89326 003

Supplement 001 dated 2015-01-13

For the product(s)/product category (ies):

Production and Distribution of
Vacuum Blood Collection Tubes with Needles,
Disposable Blood Collection Needles,
A.V. Parafix Needles, Disposable Syringes,
Auto Disposable Syringes,
Safety Syringes with Retractable Needles,
Disposable Infusion Sets,
Disposable Blood Transfusion Sets,
Scalp Ven Sets, Disposable Hypodermic Needles,
Disposable Dental Needles, Disposable IV Catheters,
Disposable Insulin Syringes,
Insulin Pen Needles, Heparin Caps,
Three-way Stopcocks (with Extension Tubes),
Blood Lancet, Sterile Surgical Blade,
Disposable Scalpel with Plastic Blade,
Tracheostomy Tubes for Single Use,
Foley Catheter (catheter), Endotracheal Tubes,
Suction Connecting Tube (with and without Yankaeur Handle),
Oxygen Masks, Nebulizer Masks,
Nebulizers, Venturi Masks, Nasal Oxygen Cannula,
Laryngeal Mask Airway, Silicone Urinary Catheter,
Latex/Oropharyngeal Airway, Sterile Isolation Catheter,
Sterile Section Catheter, Sterile Rectal Tube,
Sterile Urine Drainage Bag, Drainage Tube,
Sterile Feeding Tubes, Disposable Stomach Tubes,
Disposable Umbilical Cord Clamps,
Disposable I.D. Brackets,
Disposable Mucous Extractors,
Disposable Gynecologic Specimens,
Disposable Gynecologic Suits,
Sterile Vaginal Speculum,

Page 2 of 3

TÜV SÜD Product Service GmbH - Zertifizierungsstelle - Referenz-Nr.: 8028 München - Germany



Attachment for Certificate No. Q2M 14 09 89326 004

Supplement 001 dated 2015-01-13

Disposable Dressing Kits,
Sterile Latex Examination Gloves,
Sterile Latex Surgical Gloves,
Gauze Swab (Sponges),
Gauze Bandages, Lap Sponges,
Non-woven Sponges (with X-ray),
Sterile Parafix Gauze Dressing, Gauze Balls,
Absorbent Gauze (Gauze in Roll, Gauze in Zigzag, Cutted Gauze),
Absorbent Cotton Balls, Absorbent Cotton Balls,
Absorbent Zigzag Cotton, Eye Pad,
Alcohol Swabs / Pads, First Aid Kits,
Bandages, Elastic Bandages, Crepe Cotton Bandages,
Fast Dressing Bandages (1st Aid Bandages),
Pouch of Parafix Bandages, Cast Paddings,
Sterile Wound Dressing, Sterile Wound Dressing,
Zinc Oxide Adhesive Plaster,
Surgical Tapes (Non-woven, PE, Silk, Fabric),
Sports Tape (Include Kinesiology Tape, Rigid Tape and Fix Tape),
Surgical Gowns, Face Masks, Shoe Covers, Goggles,
Underpads, Surgical Towels, Surgical Drapes,
Isolation Gowns, Protective Clothing,
Non-woven caps, Medical Bed Sheets,
Cotton Swab (Cotton Tipped Applicators),
Tongue Depressor,
Digital Thermometers,
Infrared Ear Thermometers,
Digital Blood Pressure Monitors (Digital Sphygmomanometers),
Aneroid Sphygmomanometers, Stethoscopes
Mueh, CRT3, 2015-01-14

Hans-Mueller Jucker

Page 2 of 3

TÜV SÜD Product Service GmbH - Zertifizierungsstelle - Referenz-Nr.: 8028 München - Germany





HEADQUARTERS:

Henso Medical (Hangzhou) Co.,Ltd.
2-818 Zhongtian MCC, Tongpu Rd.,
310012 Hangzhou, China



+86 571 8604 3296



+86 571 8604 6389

HENSO (Hong Kong):

Healthaw Medical Limited
Suite 01 10/F, Carnival Commercial Building,
18 Java Road, North Point, HK



info@healthaw.com