

KT-130A

ELITE I single head stethoscope 精选I型A

- Special chestpiece design ,zinc alloy head with chrome plated
- **Kindcare**® design diaphragm ,LOGO print available
- **Kindcare**® Patent ear tips
- Inner spring binaural and 1 piece molded PVC tubing



KT-130B

ELITE I single head stethoscope 精选I型B

- Same as KT-130A, only with matte satin finish



KT-130C

ELITE I single head stethoscope 精选I型C

- Same as KT-130A, only with matte gold plated finish

Ear tips for choice



A

B

C



KT-140A

ELITE II cardiology stethoscope 精选II 型A

- Special chestpiece design ,zinc alloy head with chrome plated
- **Kindcare**® design diaphragm , logo print available
- **Kindcare**® Patent ear tips
- Binaural tubing assembly for excellent acoustic perception



KT-140B

ELITE II cardiology stethoscope 精选II 型B

- Same as KT-140A ,only with matte satin finish



KT-140B

ELITE II cardiology stethoscope 精选II 型C

- Same as KT-140A ,only with matte gold plated finish

Ear tips for choice



听诊器系列

Stethoscope series



KT-104A

ELITE III adult single head stethoscope 精选III型A

- Special chestpiece design , zinc alloy head with chrome plated
- **Kindcare** design diaphragm , LOGO print available
- Inner spring binaural and 1 piece molded PVC tubing

KT-104B

ELITE I single head stethoscope 精选III型B

- Same as KT-104A , only with matte satin finish



KT-104C

ELITE I single head stethoscope 精选III型C

- Same as KT-104A , only with matte gold plated finish



**KT-107****Stainless steel adult stethoscope(不锈钢成人听诊器)**

- For adult, deluxe type, with inner spring binaural tubing assembly,
- Special binaural design to suit ear structure
- Diameter of head : 47.5mm
- Foam box packing

KT-108**Stainless steel child stethoscope(不锈钢儿童听诊器)**

- For child, deluxe type, with inner spring binaural tubing assembly,
- Special binaural design to suit ear structure
- Foam box packing

**KT-109****Cardiology stainless steel stethoscope(不锈钢心脏科听诊器)**

- For adult, cardiology type, feature a deep cone-shaped bell
- Chestpiece comes with non-chill ring
- With inner spring binaural tubing assembly
- Special binaural design to suit ear structure
- With uniquely designed binaural thick tubing assembly
- Foam box packing

Ear tips for choice

A

B

C

**KT-118****Cardiology master stethoscope(精确型听诊器)**

- Deluxe zinc single head, special designed head for better sound correct,
- Light copper binaural, Chestpiece comes with non-chill ring, **Kindcare** patented soft ear-tips
- With inner spring binaural tubing assembly,
- Special binaural design to suit ear structure
- With uniquely designed binaural thick tubing assembly
- Foam box packing

KT-120**Stainless steel CLASS II stethoscope (成人不锈钢听诊器)**

- Deluxe stainless steel head, Light copper binaural
- Chestpiece comes with non-chill ring, **Kindcare** patented soft Ear-tips
- With inner spring binaural tubing assembly
- Special binaural design to suit ear structure
- Foam box packing

**KT-119****Cardiology Class III stethoscope(多功能不锈钢听诊器)**

- Deluxe stainless steel head, one side for child, the other side for adult
- Light copper binaural, Chestpiece comes with non-chill ring, **Kindcare** patented soft ear-tips
- With inner spring binaural tubing assembly
- Special binaural design to suit ear structure
- With uniquely designed binaural thick tubing assembly
- Foam box packing





KT-100A


Single head stethoscope for Adult(成人单头听诊器)

- Silver head (Oxidized & Anodized head for choice)
- Standard Y-tubing, Diameter of head :46mm

KT-100D

Single head stethoscope for Child(儿童单头听诊器)

- Diameter of head :37mm



KT-100B

Colored single head stethoscope for Adult(成人单头彩色听诊器)

- Attractive colored head with matching color tubing
- Diameter of head: 46mm

KT-100D1

Colored single head stethoscope for Child(儿童单头彩色听诊器)

- Diameter of head: 38mm



KT-100C

Single head stethoscope for Adult(成人单头听诊器带防冷圈)

- Aluminum chestpiece with a wide range of tubing colors that match non-chill rings
- Diameter of head :46mm

KT-100C1

Single head stethoscope for Child(儿童单头听诊器带防冷圈)

- Diameter of head :38mm



KT-100E

Special single head stethoscope for Adult(成人单头听诊器)

- Silver Anodized head, special design
- Standard Y-tubing, Diameter of head :46mm

**KT-101A**Dual head stethoscope for Adult(成人双头听诊器)

- Silver head (oxidized & anodized head for choice)
- Standard Y-tubing, diameter of head: 47mm

KT-101DDual head stethoscope for Child(儿童双头听诊器)

- Diameter of head: 38mm

**KT-101B**Colored dual head stethoscope for Adult(成人双头彩色听诊器)

- Attractive colored head with matching color tubing
- Diameter of head: 47mm

KT-101D1Colored dual head stethoscope for Child(儿童双头彩色听诊器)

- Diameter of head: 38mm

**KT-101C**Dual head stethoscope for Adult(成人双头听诊器带防冷圈)

- Aluminum chest piece with a wide range of tubing colors that match non-chill rings
- Diameter of head: 47mm

KT-101C1Dual head stethoscope for Child(儿童双头听诊器带防冷圈)

- Diameter of head: 38mm

听诊器系列

Stethoscope series

KT-102A

The sprague rappaport type stethoscope(多功能听诊器)

- Silver color head, standard tube
- Comes complete with a clear multi-pocket bag for relevant accessories and replacement parts
- Various colors available
- Clear PVC zipper case or box packing

KT-102C

The sprague rappaport type stethoscope(多功能听诊器)

- With special PVC protective ring



KT-102B1

The sprague rappaport type stethoscope(半透明 多功能听诊器)

- Silver color head, transparent tube
- Other fitting is same as KT-102A

KT-102B2

The sprague rappaport type stethoscope(彩色 多功能听诊器)

- Attractive colored head with matching color tubing
- Other fitting is same as KT-102A



KT-102B3

KT-102B3 Sprague rappaport type stethoscope(多功能听诊器)

- Attractive gold-plated head and binaural, more deluxe
- Other fitting is same as KT-102A



KT-103**Dual head special type for Adult(专用型双头听诊器)**

- With inner spring binaural tubing assembly
- Special binaural design to suit ear structure
- Comes with inner-spring binaural tubing
- Clear PVC zipper case or box packing

KT-105A**Dual head deluxe type for Adult(豪华型成人双头听诊器)**

- Aluminum chestpiece with a wide range of tubing colors that match non-chill rings
- Diameter of head: 47mm
- Comes with inner-spring binaural tubing
- Special binaural design to suit ear structure
- Clear PVC zipper case or box packing

KT-105B**Dual head deluxe type for Child(豪华型儿童双头听诊器)**

- Diameter of head: 38mm

**KT-111(ROUND SHAPE)****Acrylic stethoscope(亚克力听诊器)**

- Chestpiece in acrylic material with polishing finish
- Inner-spring binaural tubing
- PVC soft eartips
- PVC zipper case or box packing

**KT-501(HEART SHAPE)****Acrylic stethoscope(亚克力听诊器)**

- Heart shape chestpiece, other fitting is same as KT-501



KT-106

Stethoscope for Baby and Neonatal(新生儿听诊器)

- KT-106A Aluminum chestpiece
- KT-106B zinc chestpiece
- KT-106C stainless steel chestpiece
- For neonatal and baby use
- Tubing available in attractive colors, normal binaural

KT-114

Special single head stethoscope(专用型单头听诊器)

- Zinc head, powder-coating finish,
- Comes with non-chill ring,
- With inner spring binaural tubing assembly,
- Uniquely designed for the shape of head,
- Special binaural design to suit ear structure
- PVC zipper case or box packing



KT-110

Bowles stethoscope(可转动式听诊器)

- Zinc chestpiece with chrome-plated silver finish
- With inner spring binaural tubing assembly
- Uniquely designed for the shape of head
- Special binaural design to suit ear structure
- PVC zipper case or box packing

KT-113A

Special dual head stethoscope(多用型听诊器)

- Chrome-plated binaurals and fully rotatable chestpiece for adult and child.
- Internal spring binaural for professional purpose.
- PVC zipper case or box packing



KT-113B

Special dual head stethoscope(多用型听诊器)

- Chrome-plated binaurals and fully rotatable chestpiece for adult and child.
- Internal spring binaural for professional purpose.
- PVC zipper case or box packing

KT-122

ELITE type 3 set stainless steel stethoscope(三件套听诊器)

- With 3 type head : adult ,child ,pediatric



听诊器系列

Stethoscope series

KT-115

Stethoscope with clock(豪华型时钟听诊器)

- Aluminum chest piece with a wide range of tubing colors that match non-chill rings
- Diameter of head: 47mm
- Comes with inner - spring binaural tubing
- Special binaural design to suit ear structure
- The head with clock from JAPAN
- Clear PVC zipper case or box



KT-116

Sprague rappaport stethoscope(多功能时钟听诊器)

- Silver color head, standard tube
- Comes complete with a clear multi-pocket bag for relevant accessories and replacement parts
- Various colors available
- Clear PVC zipper case or box packing



KT-117

Special first aid stethoscope(户外听诊器)

- Special chestpiece for outdoor use
- Internal spring binaural for professional purpose
- PVC zipper case or box packing

听诊器系列

Stethoscope series

KT-112

Frequency conversion single head stethoscope(豪华单头听诊器)

- Zinc alloy chestpiece with chrome-plated silver finish
- PVC material membrane with special-designed central hole for better sound transfer
- Inner-spring binaural tubing, PVC soft eartips
- PVC zipper case or box packing



KT-121

Stainless steel single head stethoscope(不锈钢成人单头听诊器)

- For adult, deluxe single head type, with inner spring binaural tubing assembly
- Special binaural design to suit ear structure
- Diameter of head: 47.5mm
- Foam box packing



KT-401

Fetal stethoscope(胎声听诊器)

- For pregnancy use



KT-402

Fetal stethoscope(多用途胎声听诊器)

- For pregnancy use and normal adult use

听诊器系列

Stethoscope series



KT-201

Single head stethoscope for teaching use(单头教学用听诊器)

- With two binaurals, PVC zipper case or box packing



KT-202

Dual head stethoscope for teaching use(双头教学用听诊器)

- With two binaurals, PVC zipper case or box packing



KT-203

Dual head stethoscope for teaching use(双头教学用听诊器)

- With two binaurals, non-chill rings
- PVC zipper case or box packing



KT-403

Fetal stethoscope(胎声听诊器)

- For pregnancy use



KT-301

Plastic stethoscope for adult(塑料玩具听诊器)

- Full set are all produced by plastic